

LGT-2202DVT

Elevated VTE Prevention Experience



A Smart System that Improves VTE Prevention in Hospital Settings

LGT-2202DVT delivers comfortable air compression therapy in a lightweight and remarkably compact package to effectively address VTE problems of post-surgical patients. It is designed and developed to further simplify the workflow and expand clinical capabilities of VTE prevention in hospital settings.

Benefits of the treatment

- VTE prevention
- Circulation enhancement
- Acute and chronic edema reduction
- Improve the healing rate of leg ulcers
- Venous stasis treatment
- Venous insufficiency prevention and treatment





Automatic Garment Recognition



One-Button-Start System



Real-Time Monitoring



Quiet Operation

FEATURES

- Compact dimensions and lightweight
- Battery-operated
- Quiet operation
- Automated garment recognition
- Intuitive user interface
- Alternate and synchronous inflation
- Single leg treatment supported
- Vascular refill detection, safe and smart compression cycles
- Real-time pressure monitoring
- Gradient compression
- High-quality connector, easy to connect and disconnect
- Complete range of garments, made of high quality and breathable fabrics





COMPLETE RANGE OF GARMENTS

Order No.	Sleeve	Size	Proximal Circumference	
220094	Calf Sleeve	S	<40.6cm	 <p>Calf Sleeve</p>
220094M		M	<48.3cm	
220094L		L	<55.9cm	
220095M	Foot Sleeve	M	<40.6cm	 <p>Foot Sleeve</p>
220092	Thigh Length Sleeve	S	<55.9cm	 <p>Thigh Length Sleeve</p>
220092M		M	<71.1cm	
220092L		L	<91.4cm	
220093M	Calf and Foot Sleeve	M	<53.3cm	 <p>Calf and Foot Sleeve</p>

Health Life · Longest Care

Longest
PHYSIO & REHAB DEVICES

Guangzhou Longest Medical Technology Co., Ltd.

Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-020-66353999

Fax: 86-020-66353920

E-mail: export@longest.cn

Website: www.longestmedical.com

