Certification



Technical Specification

Pressure Range: 1.0~4.0 Bar (Step by 0.1 Bar)
Frequency: 1~17Hz
Pre-set Protocol: Yes
Circumference: Yes
Dimension: 39.5 (L) x 39 (W) x 19.4 (H) cm
Standard Weight: 11.3kg
Langages: English, Italian, German, Spanish, French

Longest Guangzhou Longest Medical Technology Co., Ltd.

Address: 301&401 of Building 2& Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-20-66353999 Fax: 86-20-66353920


E-mail: export@longest.cn Website: www.longestmedical.com



PowerWave LGT-2500S

The ideal body shaping treatment Acoustic waves provide an ideal body shaping treatment option











• Applications

- ► Anti-cellulite treatment
- Body reshaping
- ► Skin elasticity improvement
- ► Treatment of skin irregularity after liposuction
- ► Connective tissue tightening
- ► Scar and wrinkle smoothing
- ► Treatment of stretch marks
- ► Lymphatic system stimulation














• Effects

- ► Body contouring



PowerWave LGT-2500S

Acoustic Wave Therapy Device

Acoustic wave treatment uses acoustic sound waves together with the body's own ability to repair tissue by stimulating circulation in all layers of the skin, including blood vessels, connective tissue and muscles.

AWT tightens and forms new connective tissue, reducing cellulite and improving skin appearance.

Benefits

- ► The acoustic wave is able to massage the tissues in depth with considerable beneficial effects.
- ► Stimulation of blood circulation and lymphatic
- ► Raising the heat tissue and activation of the elimination of fat cells
- ► Activation of metabolism

► Improvement in skin elasticity

Stimulating the production of collagen, improv-ing elasticity and texture of the skin