


LGT-2410S

Localized Cryotherapy System



EXPAND RECOVERY POTENTIAL

LGT-2410S combines the effects of cold and pressure. In the cylinder, the CO₂ liquid is pressurized. When treating, the valve is open and the gas streams out through the handpiece and hits the treated area at a temperature of -78°C and a pressure of 50 bars.

At room temperature, hitting the treated area at a distance of 4cm can lower the skin temperature of the treated area to 4 °C in the 30s.



*The cylinder is not provided.
Please purchase locally.



FEATURES

LGT-2410S is a cutting-edge localized cryotherapy device featuring chic appearance, robust structure, easy operation, remarkable durability, and superior performance. It is designed to simplify workflows, improve productivity, and create satisfying experience.

Intuitive Touchscreen

- Clear and comfortable view, easy control
- Therapeutic encyclopedia with image guidance
- Real-time skin temperature and remaining treatment time display

Rotary Knob

- Quick setting of parameters

Blue Light

- Positioning the treated area for more precise treatment
- Kill bacteria and promote healing

Monitoring Screen

- Give real-time feedback of skin temperature and remaining treatment time for more controlled treatment

Ergonomic Handpiece


- Comfortable to use
- Anti-frosting design of the nozzle for long service life

Integrated Trolley

- Greater maneuverability
- Large storage space



EXTENSIVE APPLICATIONS



Benefits of the treatment

- Immediate pain relief
- Reduced inflammation
- Accelerated recovery
- Muscle relaxation
- Improved sports performance
- Skin rejuvenation



Health Life · Longest Care

Longest
PHYSIO & REHAB DEVICES

Guangzhou Longest Medical Technology Co., Ltd.

Address: 301 & 401 of Building 2 & Building 3, No.96, Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, Guangdong Province, 511399, China

Tel: 86-020-66353999

Fax: 86-020-66353920

E-mail: export@longest.cn

Website: www.longestmedical.com

